



Z I V I D

Accurate and high-speed 3D vision
for collaborative robotics

Industry 4.0 robotics and 3D vision

Symbiotic and synergistic relationship



Advancements in 3D vision drives increased productivity and growth

Speed and reliability for existing applications

Enabling new applications

Ease of deployment, operation and flexibility





Bigger robots
- bigger objects

bigger objects
- less accurate vision

Smaller robots
- smaller objects

smaller objects
- more accurate vision





Smaller robots - lots of smaller parts

The need for highly accurate 3D
vision in collaborative robotics

detecting | picking | placing

Question #1

How accurate 3D vision
do I need?

Question #2

Which of these 3D cameras
meet my needs?





Question #2

Which of these 3D cameras
meet my needs?

3D cameras are measurement instruments
We need to step up
our game...



ISO 5725

Common terminology on accuracy

Precision


Describing **random errors**,
a measure of statistical variability.

Trueness

Describing **systematic errors**,
a measure of statistical bias.

Accuracy


Describing a the combination of
random and systematic errors.
Sum of Precision and Trueness.




ISO 5725 and 3D cameras

Precision, trueness and accuracy applied

Lower accuracy



Higher accuracy



Common terminology on metrics


Accuracy on what?




Point



Plane



Dimension



Sphere

Industry standard technical specifications

Reproducible specifications

PRODUCT DESCRIPTION

LOCTITE® 270™ provides the following product characteristics.

Technology	Anaerobic
Chemical Type	Dimethyl acrylate ester
Appearance (uncured)	Green liquid
Fluorescence	Fluoresces under UV light
Applications	Steel, steel/steel - requires no mixing
Viscosity	Low
Cure	Anaerobic
Secondary Cure	Activator
Application	Threadlocking
Strength	High

This Technical Data Sheet is valid for LOCTITE® 270™ manufactured from the dates outlined in the "Manufacturing Date Reference" section.

LOCTITE® 270™ is designed for the permanent locking and sealing of threaded fasteners. The product cures when confined in the absence of air between close fitting metal surfaces and prevents loosening and leakage from shock and vibration. LOCTITE® 270™ is particularly suited for heavy duty applications such as studs into motor housings, nuts onto studs in pump housings and other fasteners where high strength is required. LOCTITE® 270™ provides robust curing performance. It not only works on active metals (e.g. brass, copper) but also on passive substrates such as stainless steel and plated surfaces. The product offers high temperature performance and oil tolerance. It tolerates minor surface contamination from various oils, such as cutting, lubrication, anti-corrosion and protection fluids.

NSF International Registered to NSF Category P1 for use as a sealant where there is no possibility of food contact in and around food processing areas. **Note:** This is a regional approval. Please contact your local Technical Service Center for more information and clarification.


TYPICAL PROPERTIES OF UNCURED MATERIAL

Specific Gravity @ 25 °C 1.1
 Flash Point 500MSDS
 Viscosity (Brookfield - RVT, 25 °C, mPa s (cP))
 Spindle 2 speed 20 rpm 400 to 600
 Viscosity (Cone & Plate, 25 °C, mPa s (cP))
 Cone 060.1 Ti @ shear rate 129 s 450

TYPICAL CURING PERFORMANCE


Cure Speed vs. Substrate

The rate of cure will depend on the substrate used. The graph below shows the breakaway strength developed with time on M10 steel nuts and bolts compared to different materials and tested according to ISO 10964.



Cure Speed vs. Bond Gap

The rate of cure will depend on the bondline gap. Gaps in threaded fasteners depends on thread type, quality and size. The following graph shows shear strength developed with time on steel pins and collars at different controlled gaps and tested according to ISO 10123.



Conditions


- Camera settings
- Working distance
- Ambient light
- Ambient temperature

Technique

- Test-setup
- Test procedure
- Calculations of result

Point metrics

Zivid One+ M - typical specifications



Point precision

$L_a = 750$ lux

$D = 1600$ mm

$D = 1600$ mm, $P_b = 1.8$ x

60 μ m

70 μ m

570 μ m

480 μ m

Conditions, unless otherwise specified:

Ambient temperature (T_a) = 25C

Ambient light (L_a) = 200 lux


Aperture (A) = f/5.6

Projector brightness (P_b) = 1x


Imager gain (I_g) = 1x

Point metrics

Zivid One+ M - typical specifications



Point precision vs. working distance



Conditions, unless otherwise specified:

Ambient temperature (T_a) = 25C

Ambient light (L_a) = 200 lux


Aperture (A) = f/5.6

Projector brightness (P_b) = 1x

Imager gain (I_g) = 1x

Global planarity metrics

Zivid One+ M - typical specifications



Global planarity trueness	90 um
D = 1600mm	670 um
D = 1600mm, Pb = 1.8x	690 um

Global planarity accuracy	140 um
D = 1600mm	1030 um
D = 1600mm, Pb = 1.8x	1000 um

Conditions, unless otherwise specified:

Ambient temperature (Ta) = 25C

Ambient light (La) = 200 lux

Aperture (A) = f/5.6

Projector brightness (Pb) = 1x

Imager gain (Ig) = 1x



Need for highly accurate 3D vision in collaborative robotics

Small robots – many small object – highly accurate vision for detection, picking and placing

- Understanding accuracy of a camera is tricky
- Lack of common terminology, conditions, set-up and technique
- 3D cameras are measurement instruments

We need to step up our game on specifying accuracy

Enabling a system developer assess whether a 3D camera meets the accuracy requirements of the application

- ISO 5725 framework for Precision, Trueness and Accuracy
- Point, plane, dimension and sphere metrics
- Reproducible specifications with industry standard specifications with conditions, set-up and technique

Next steps

Zivid One+ datasheet

zivid.com/downloads

See for yourself!

- Small
- Medium
- Large

zivid.com/zivid-one-plus-small-3d-camera

zivid.com/zivid-one-plus-medium-3d-camera

zivid.com/zivid-one-plus-large-3d-camera

Free 1:1 online demo

zivid.com/schedule-a-free-zivid-demo

ORDER

Learn more at
www.zivid.com

Human-like **vision** for robots

